





icare
HOME

TA022-088-1.0

EN	INSTRUCTION MANUAL FOR HEALTH CARE PROFESSIONALS
DE	BEDIENUNGSANLEITUNG FÜR MEDIZINISCHES FACHPERSONAL
ES	MANUAL DE INSTRUCCIONES PARA PROFESIONALES SANITARIOS
IT	MANUALE DI ISTRUZIONI PER PROFESSIONISTI SANITARI
SV	ANVÄNDARHANDBOK FÖR VÅRDPERSONAL
FI	KÄYTTÖOHJE TERVEYDENHUOLLON AMMATTILAISILLE
NO	BRUKSANVISNING FOR HELSEPERSONELL
PL	INSTRUKCJA OBSŁUGI DLA PRACOWNIKÓW OPIEKI ZDROWOTNEJ
DA	BRUGERMANUAL FOR SUNDHEDSPERSONALE
PT	MANUAL DE INSTRUÇÕES PARA PROFISSIONAIS DA SAÚDE
NL	GEbruIKSAANWIJZING VOOR MEDISCH PERSONEEL
FR	MODE D'EMPLOI À L'ATTENTION DES PROFESSIONNELS DES SOINS DE SANTÉ


19


20

21

22

23

24

25

26

27

28

29

30

31

32

33

34


35


36


37


38


39


40


41

42

43

44

45

46

ENGLISH

1.	Safety instructions.....	1
2.	Indications for use.....	2
3.	Contraindications	2
4.	Warning	2
5.	Precaution	3
6.	Introduction	3
7.	Package contents	3
8.	Before you start.....	3
9.	Setting up the tonometer.....	4
	9.1 Installing or changing the battery	4
	9.2 Turning the tonometer on	4
	9.3 Loading the probe.....	4
10.	Using the tonometer.....	5
	10.1 Choosing the measurement mode	5
	10.2 Adjusting the measurement position	5
	10.3 Automatic eye recognition.....	5
	10.4 Taking the measurements	5
11.	Reading the measurement data.....	6
12.	Troubleshooting	6
13.	Replacing the probe base.....	7
14.	Cleaning the probe base	8
15.	Cleaning and disinfection.....	8
16.	Accessories.....	8
17.	Lifetime.....	8
18.	Technical and performance data	8
19.	Symbols	9
20.	Electromagnetic declaration	10
21.	Training.....	12

ENGLISH

Icare® HOME (Model: TA022) INSTRUCTION MANUAL

The information in this document is subject to change without prior notice.
In a conflict situation the English version prevails.

This device complies with: Medical Device Directive 93/42/EEC
Canadian Medical Device Regulations
RoHS Directive 2011/65/EU

Copyright © 2017 Icare Finland Oy
Made in Finland

Icare Finland Oy
Äyritie 22, 01510 Vantaa
Tel. +358 9 8775 1150, Fax +358 9 728 6670

0598

www.icaretonometer.com, info@icarefinland.com

1. SAFETY INSTRUCTIONS

WARNING!

Do not push the tonometer into the eye (the tip of the probe should be 4-8mm, or 5/32-5/16", from the eye).

WARNING!

Keep the tonometer out of the reach of children, because the probe base, battery compartment cover and probes are so small that a child could swallow them.

WARNING!

The probe tips are for single-use only, and are packaged sterile.

WARNING!

To prevent contamination, do not touch the bare probe, do not use a probe if it touches a non-sterile surface like a table or a floor.

WARNING!

Health care professionals must inform patients not to modify or discontinue their treatment plan without receiving instructions from the health care professional.

WARNING!

For cybersecurity do not connect to the USB port except when uploading patient measurement data. Also the tonometer does not allow you to take any measurements when the USB is connected. All other traditional cybersecurity controls (anti-virus software, malware software, separate network for the device, etc.) do not apply since the device is stand alone, is not networked, and does not contain operating system software.

WARNING!

Do not change the batteries or probe base when the USB cable is connected.

WARNING!

No modification of this equipment is allowed.

WARNING!

Use only the original and certified probes made by the manufacturer. The probes are for single-use (single pair of measurement sequences) only. Use probes taken only from the intact, original packaging. The manufacturer cannot guarantee sterility of the probe once the seal is compromised. Re-sterilization or re-use of the probe could result in incorrect measurement values, in the breakdown of the probe, cross-contamination of bacteria or viruses, and infection of the eye. Re-sterilization or re-use will void all responsibilities and liabilities of the manufacturer concerning the safety and effectiveness of the tonometer.

WARNING!

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

**WARNING!**

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

**PRECAUTION!**

- *When you have opened the package, check for any external damage or faults, particularly for damage to the case. If you suspect that there is something wrong with the tonometer, contact the dealer who sold the tonometer to you.*
- *Use the tonometer only for measuring intraocular pressure. Any other use is improper and the manufacturer is not liable for any damage arising from improper use, or for the consequences of such use.*
- *Never open the casing of the tonometer, except for the battery compartment.*
- *Never allow the tonometer to get wet.*
- *Do not use the tonometer near flammable substances, including flammable anesthetic agents.*
- *Certain microbiological agents (for example, bacteria) can be transmitted from the forehead or cheek support. To prevent this, clean the forehead and cheek support for each new patient with disinfectant. See the chapter 'Cleaning and disinfection'.*
- *The tonometer conforms to EMC requirements (IEC 60601-1-2), but interference may occur within the tonometer if used near (<1m) a device causing high intensity electromagnetic emissions, such as a cellular phone. Although the tonometer's own electromagnetic emissions are well below the levels permitted by the relevant standards, they may cause interference in other, nearby devices, for example sensitive sensors.*
- *If you do not use the tonometer for a long time, remove the batteries, as they may leak.*
- *Be sure to dispose of the single-use probes properly (for example, in a bin for metal waste).*
- *Batteries, packaging materials and probe bases must be disposed of according to local regulations.*
- *Make sure you use batteries with built-in PTC protection, for example Energizer Lithium Photo 123 3V CR123A.*
- *Do not cover the eye recognition transmitters or sensor during the measurement, for example with your fingers. Keep your hand, hair etc. and objects such as pillows away from the temple side of your eye, as they produce an infrared reflection that causes an error.*
- *The tonometer turns off automatically after 3 minutes if you do not use it.*
- *Update the tonometer's time to your local time. It is done automatically by performing the steps 1 and 2 under section 11. Reading the measurement data.*

2. INDICATIONS FOR USE

The Icare HOME tonometer is a prescription device intended for monitoring of intraocular pressure (IOP) of the human eye. It is indicated for use by patients or their caregivers under supervision of an eye care professional.

3. CONTRAINDICATIONS

Patients with any of the following conditions should not be prescribed the HOME tonometer:

1. Active ocular infection including infectious conjunctivitis
2. Recent ocular trauma including corneal laceration or corneal/scleral perforation
3. Disabling arthritis or limited motor coordination affecting self-handling of the Icare tonometer
4. Blepharospasm
5. Nystagmus

4. WARNING

Patients with the following conditions are generally not eligible for use of the HOME tonometer as they have the potential to introduce unsafe conditions during use or to impair measurement acquisition:

1. Uncorrected near visual acuity of 20/200 or worse
2. Only one functional eye
3. Poor or eccentric fixation
4. Hearing impairment to the extent that the individual cannot hear and converse with others without an assistive aid and/or sign language

5. Contact lens use
6. Dry eyes
7. Keratoconus
8. Microphthalmos
9. Buphthalmos
10. Cataract extraction within last 2 months

5. PRECAUTION

The safety and effectiveness of the Icare HOME tonometer has not been evaluated for patients with:

1. High corneal astigmatism >3d
2. History of prior incisional glaucoma surgery or corneal surgery including corneal laser surgery
3. Corneal scarring
4. Central corneal thickness greater than 0.60mm or less than 0.50mm
5. Known history of difficulty in obtaining Goldmann IOP measurements or any factors that might contribute to inaccurate Goldmann IOP measurements (e.g. lid squeezing or tremor)

6. INTRODUCTION

The Icare HOME tonometer is a hand-held device for self-use. The great advantage is that a topical anesthetic is not needed.

The tonometer uses the rebound method. A small and light single-use probe makes contact with the eye very briefly. The tonometer measures the deceleration of the probe and the rebound time, and calculates the IOP from these parameters.

A measurement sequence includes six measurements. The probe moves to the cornea and back during every measurement. As a result, after the six measurements the tonometer calculates the final IOP and stores it with other information in the tonometer's memory, including date, time, eye identification (right or left) and measurement quality.

The Icare HOME tonometer can record over one thousand measurement results. You can copy the recorded measurement information to a PC through a USB cable for management of your glaucoma patients.

7. PACKAGE CONTENTS

The package contains:

- Icare HOME tonometer
- 10 sterilized single-use probes
- 2 batteries
- USB memory stick including the instruction manual for health care professionals and the Icare LINK software
- USB cable for connecting the Icare HOME tonometer to a PC with Icare LINK software
- Instructions for downloading the Icare LINK software and registration of the tonometer
- Patient guide
- Support position tags
- Warranty card
- Carrying case
- Wrist strap
- Probe base cleaning container

8. BEFORE YOU START

Find the main parts, buttons and indicator lights of the tonometer. The figures are located at the beginning of this document.

FRONT PARTS (Figure 1)

6. Probe base incorporating indicator light
7. Eye recognition transmitter
8. Eye recognition sensor
9. Cheek support
10. Forehead support

TOP PARTS (Figure 2)

11. Measurement button
12. Forehead support position indicator

SIDE PARTS (Figure 3)

13. Forehead support adjustment wheel
14. Cheek support adjustment wheel

BOTTOM PARTS (Figure 4)

15. Battery cover
16. Silicone lid (USB cover)
17. Type label

BACK PANEL (Figure 5)

18. Measurement button
19. LOAD light
20. MEASURE light
21. REPEAT light
22. DONE light
23. SERVICE light
24. BATTERY light
25. POWER button

9. SETTING UP THE TONOMETER

Setting up your Icare HOME tonometer is easy, with few steps. The following subchapters describe how you get started.

9.1 INSTALLING OR CHANGING THE BATTERY

Lift the silicone lid that protects the USB port and keeps the battery compartment cover in place. Open the battery compartment cover by pressing the silicone lid slightly and sliding the battery compartment cover as shown in the Figure 6.

1. Silicone lid
2. Battery cover

Insert two CR123A lithium batteries in the correct order: (+) end upwards as shown in the figure left. Close the cover firmly and press the silicone lid in place to cover the USB port.

9.2 TURNING THE TONOMETER ON

Press the power button (20) to turn the tonometer on. The lights (14-19) are turned on briefly. Following a brief pause, the Load light flashes on the back panel to remind the user to load the single-use probe into the tonometer prior to measurement.

9.3 LOADING THE PROBE

The Icare HOME tonometer uses single-use probes that are packed in a plastic tube and wrapped in blister packs as shown in the Figures 8 and 9.

To load the probe:

1. Unwrap the probe (Figure 10).
2. Remove the lid of the probe container as shown in the Figure 11. Point the tonometer upward.
3. Drop the probe into the probe base by turning the probe container upside down as shown in the Figure 12.
4. Press the measurement button briefly (1 second) to activate the probe (Figure 13).
5. The probe moves rapidly back and forth (Figure 14).
6. If the green Measure light flashes, the probe is loaded correctly and ready for measurement (Figure 15).

10. USING THE TONOMETER

10.1 CHOOSING THE MEASUREMENT MODE

The tonometer can operate in two modes:

Series mode

The series mode is especially useful in self-tonometry. In the series mode, keeping the button pressed down (see Figure 16) initiates the measurement function, and the tonometer takes six rapid measurements one after the other to obtain the final IOP reading. The button must be pressed for at least 3 seconds.

Single mode

You can use the single mode to take individual measurements one at a time. The single mode is especially useful for those patients who tend to blink heavily. Here you press the measurement button briefly (1 second) for each of the six measurements to obtain the final IOP reading (Figure 17).

10.2 ADJUSTING THE MEASUREMENT POSITION

The tonometer has two adjustable supports (see points 4 and 5 in Figure 1), one for the forehead and one for the cheek, as shown in Figure 18. The supports are for ensuring accurate measurement distance and alignment.

To adjust the measurement position for your patient:

1. Adjust the supports using the adjustment wheels as shown in Figure 18.
2. Keep the probe horizontal and pointing perpendicularly to the center of the cornea.
3. Set the distance between the tip of the probe and the center of the cornea to be 4-8mm (5/32-5/16") as shown in Figure 19.
4. Read the distance setting (forehead A•1, A•2, etc., cheek B•1, B•2, etc.) between the arrows on the scale (see point 7 in Figure 2) of the supports (see Figure 20), and write it down on a support position tag for the patient.
5. Do the same for the other eye as well unless only one eye needs monitoring.
6. Verify that the support positions are correct each and every time the patient comes in for a clinic visit.

10.3 AUTOMATIC EYE RECOGNITION

The tonometer includes an automatic eye recognition system that identifies which eye, right or left, you are measuring. The system has two infrared LED transmitters just below the probe base and one infrared LED sensor above the probe base, as in Figure 21. The right-hand transmitter sends invisible infrared light to the right and the left-hand transmitter to the left. The infrared light reflects from your nose to the sensor. The sensor knows from which transmitter the reflected infrared light came, and thus which eye you are measuring. The resulting eye indication is included in the data that you can transfer to a PC, as described in the section 11.

EYE RECOGNITION COMPONENTS (Figure 21)

1. Left and right infrared transmitter.
2. Infrared sensor.

10.4 TAKING THE MEASUREMENTS

The probe will make a gentle and brief contact with the eye when you take the measurement. No topical anesthetic is needed. The recommended frequency of measurements is 3-4 daily with a maximum of 5-6.

To measure intraocular pressure:

1. Check that the **Measure** light still flashes on the back panel.
2. If the **Measure** light does not flash, press the **power button** and wait until the Measure light illuminates again.
3. The patient should look straight ahead at a specific point while keeping eyes wide open as shown in Figure 22.

4. Bring the tonometer near the eye, the probe pointing perpendicular to the center of the cornea without a vertical or horizontal tilt. The position is correct when the probe base light is green and appears symmetrically in the center of the patient's view. See Figure 23 and 24.
 - Correct alignment of the tonometer (see Figure 25).
 - Incorrect alignment of the tonometer (see Figure 26).
 - Incorrect alignment of the tonometer. Readjust so that you see only the front of the tonometer and the green light symmetrically in the center of your view (see Figure 27).
5. Press the measurement button:

Single mode:

Press the button briefly (1 second) and you hear a short beep, repeat it to take one measurement at a time till you hear a long beep and see the Done light (17) illuminated on the back panel (Figure 28).

Series mode:

Keep the measurement button down to obtain the sequence of six measurements till you hear a long beep and see the Done light illuminated on the back panel. The button must be pressed for at least 3 seconds (Figure 29).

6. If both eyes are measured repeat steps 1-5 using your other eye.
7. If an error occurs, press the Measure button briefly (1 second) and continue the measurement. See also section 12 Trouble shooting.
8. Press the power button for three seconds to turn the tonometer off. The probe slides out of the probe base.
9. Remove and discard the used probe properly.

11. READING THE MEASUREMENT DATA

The tonometer stores information on every complete measurement sequence of six measurements. The stored information includes the calculated final eye pressure reading in mmHg, time and date of the measurement, identification of the eye (right or left) and the quality level of the measurement. Uploading is easy:

1. Start Icare LINK software in your PC (Figure 30).
2. Connect the tonometer to the PC using the USB cable. The Load and Measure lights will flash. If no lights flash or the Service and Battery lights flash, reconnect the USB cable.
3. The internal clock of the tonometer is automatically updated to the PC's time by the Icare LINK software at this point (Figure 31).
4. Copy the data to a selected patient in the Icare LINK software (Figure 32).

More information about **Icare LINK software**
<http://www.icaretonometer.com/products/icare-link/>

12. TROUBLESHOOTING

The tonometer automatically monitors and controls the measurement position and speed of the probe during the measurements, and indicates errors with sounds and lights. The following table instructs you in error situations and explains what the different lights and sounds mean. The indicator lights are also presented in the figure below the table.

Error light	Error sound	Reason	Action
Battery (Figure 33).	No.	Battery is soon empty.	Prepare to change batteries.
Battery light is flashing.	No.	Battery is empty.	Change batteries.
Probe base light is solid red (Figure 34).	No.	Too much vertical tilt.	Press the measurement button again to clear the error message. Position the tonometer horizontally so that the probe base light is green.

Error light	Error sound	Reason	Action
Probe base light is flashing red and Measure light turns off (Figure 35).	Two long beeps.	a) Probe is too far from or too near the eye. b) Probe movement was not perpendicular to the cornea.	Press the measurement button again to clear the error message. a) Set the correct distance 4-8mm (5/32-5/16") between the probe tip and the center of the cornea. b) Set the probe perpendicular to the center of the cornea.
Repeat light is flashing and Probe base light is flashing red (Figure 36).	Two long error beeps.	a) Too much IOP deviation during the measurement, because the user did not keep the tonometer stable. b) Eye was not recognized.	Press the measurement button again to clear the error message. a) Repeat the measurement. b) Do not move the tonometer during the measurements, remove your hand or fingers from the infrared transmitters and sensor, move the patient's hair away from his/her temple side of the eye.
Service light is flashing and Probe base is flashing red (Figure 37).	Two long beeps.	Incorrect or dirty probe or probe base.	Change the probe, clean or change the probe base or contact the seller to arrange sending the device for service.

OFF (Figure 38)

1. Lights off

INITIALIZE (Figure 39)

2. LOAD (green icon)
3. MEASURE (green icon)
4. REPEAT (yellow icon)
5. DONE (green icon)
6. SERVICE (red icon)
7. BATTERY (red icon)

LOAD PROBE (Figure 40)

8. LOAD (green icon)

MEASURE (Figure 41)

9. MEASURE (green icon)

REPEAT (Figure 42)

10. REPEAT (yellow icon)

DONE (Figure 43)

11. DONE (green icon)

SERVICE (Figure 44)

12. SERVICE (red icon)

LOW BATTERY (Figure 45)

13. BATTERY (red icon)

EMPTY BATTERY (Figure 46)

14. BATTERY (red icon)

13. REPLACING THE PROBE BASE

Replace the probe base every twelve months. Change the probe, replace or clean the probe base if the Service light flashes.

Instructions for replacing the probe base:

- Turn off the tonometer.
- Unscrew the probe base collar and put it in a safe place.
- Remove the probe base by tilting the tonometer downwards and use your fingers to pull the probe base out of the tonometer.
- Insert a new probe base into the tonometer.
- Screw the collar in, to lock the probe base.

14. CLEANING THE PROBE BASE

You can reuse the probe base after careful cleaning. Clean the probe base every six months. Change the probe, clean or replace the probe base if the Service light flashes.

Instructions for cleaning the probe base:

- Fill the probe base cleaning container or other clean container with 70-100% isopropyl alcohol.
- Turn the power off.
- Unscrew the probe base collar.
- Invert the probe base over the container, drop in the probe base into the container and let soak for 5-30 minutes.
- Remove the probe base from alcohol.
- Dry the probe base by blowing clean canned or compressed air into the hole in the probe base. This will additionally remove possible residual dirt.
- Insert the probe base into the tonometer.
- Screw the collar in, to lock the probe base.

15. CLEANING AND DISINFECTION

You must clean the forehead and cheek supports for each new patient. Use a wipe dampened with a 70-100% isopropyl alcohol solution. Do not immerse the tonometer in water or other liquid. The tonometer must not be immersed or cleaned using too much water.

16. ACCESSORIES

Part number	Product Description	Weight	Dimensions
540	Probe base	4g	7 mm x 38 mm
TA022-001	Probe base collar	2g	20 mm x 15 mm
560	Wrist strap	3g	270 mm x 10 mm x 10 mm
TA022-044	Carrying case	210g	270 mm x 135 mm x 60 mm
7179	Battery cover	3g	26 mm x 23 mm x 7 mm
TA022-035	Patient guide	33g	210 mm x 90 mm x 2 mm
571	Battery 3 V, CR123A	17g	17 mm x 35 mm
TA022-037	Support position tags	40g	70 mm x 41 mm x 13 mm
575	USB cable	23g	1m
113	Probe box	55,14g	8,2 x 19,5 x 3,5 cm
543	Probe base cleaning container	3g	5,6 cm x 2 cm

17. LIFETIME

The expected service life of the device is 5 years. The maintenance described in the chapters 12-15 is required during the the expected service life. The shelf life of the probes in their intact original packaging is 3 years. We recommend that you inspect the device for mechanical and functional damage and the safety labels for legibility annually/every 12 months.

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18. TECHNICAL AND PERFORMANCE DATA

Type: TA022

Dimensions: approximately 11cm x 8cm x 3cm.

Weight: approximately 150g.

Power supply: 2 x CR123 non-rechargeable batteries (make sure you use batteries with built-in PTC protection, for example Energizer Lithium Photo 123 3V CR123A).

Measurement range: 5-50 mmHg.

Accuracy (95% tolerance interval relatively to manometry): $\pm 1,2$ mmHg (≤ 20 mmHg) and ± 2.2 mmHg (> 20 mmHg).

Repeatability (coefficient of variation): $< 8\%$.

The serial number is located on the inside of the battery compartment cover.
The lot number of the probes is on the side of the probe box and the blister packing.

There are no electrical connections from the tonometer to the patient.

The tonometer has BF-type electric shock protection.

Operation environment:

Temperature: $+10\text{ }^{\circ}\text{C}$ to $+35\text{ }^{\circ}\text{C}$

Relative humidity: 30% to 90%

Atmospheric pressure: 800hPa – 1060hPa

Storage environment:

Temperature: $-10\text{ }^{\circ}\text{C}$ to $+55\text{ }^{\circ}\text{C}$

Relative humidity: 10% to 95%

Atmospheric pressure: 700hPa – 1060hPa

Transport environment:

Temperature: $-40\text{ }^{\circ}\text{C}$ to $+70\text{ }^{\circ}\text{C}$

Relative humidity: 10% to 95%

Atmospheric pressure: 500hPa – 1060hPa

Environmental restrictions for professional use include:

- Medivac vehicles or similar where vibration or noise levels are so high that the user cannot hear error signals.

Environmental restrictions for lay operators (patients):

- Environments where noise is so high that the user cannot hear the error signals.

Mode of operation: continuous

19. SYMBOLS

Caution

See operating instructions for more information

BF-type device

Single-use disposable

Serial number

Use by <date>

Manufacturer

IP22

Protected against insertion of fingers and will not be damaged or become unsafe during a specified test in which it is exposed to vertically or nearly vertically dripping water.

Class 1 LED product

This product meets the power requirements for a Class 1 LED product to IEC/EN 60825-1 (2001) under normal operating conditions and those of single fault failure.

Keep dry

Manufacturing date

Lot number


Sterilized using irradiation

Stand by

Do not discard this product with other household-type waste. Send to appropriate facility for recovery and recycling. EU WEEE (European Union Directive for Waste of Electronic and Electrical Equipment)

20. ELECTROMAGNETIC DECLARATION

Icare HOME is class B equipment and needs special precautions regarding EMC and needs to be installed and put into service according to EMC information provided below.

Guidance and manufacturer's declaration - Electromagnetic emissions		
Icare HOME (TA022) is intended for use and should be used in the electromagnetic environment specified below.		
RF emissions CISPR 11	Group 1	Icare HOME (TA022) is battery operated and uses RF energy only for its internal function. Therefore, its RF emissions are low and are not likely to cause any interference in nearby equipment.
RF emissions CISPR 11	Class B	Icare HOME (TA022) is suitable for use in all establishments, including domestic establishments and those directly connected to public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC 61000-3-2	NOT APPLICABLE	Power level of Icare HOME (TA022) is below standard requirement, batteries are not rechargeable
Voltage fluctuations flickering emissions IEC 61000-3-3	NOT APPLICABLE	Icare HOME (TA022) batteries are not rechargeable

Guidance and manufacturer's declaration-Electromagnetic immunity			
Icare HOME (TA022) is intended for use and should be used in the electromagnetic environment specified below.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment-Guidance
Electrostatic discharge (ESD)IEC 61000-4-2	± 8 kV contact ± 2kV, ±4kV, ±8kV, ±15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast Transients/burst IEC 61000-4-4	± 2 kV 100 kHz repetition frequency	NOT APPLICABLE	Icare HOME (TA022) tonometer is not operational when connected to an external computer, Icare HOME (TA022) batteries are not rechargeable
Surge IEC 61000-4-5	±1 kV for line(s) to line(s) ±2 kV for line(s) to earth	NOT APPLICABLE	Icare HOME (TA022) tonometer is not operational when connected to an external computer, Icare HOME (TA022) batteries are not rechargeable
Voltage dips, short interruption and voltage variations on power supply lines IEC 61000-4-11	0 % UT for 0.5 cycle (1 phase) 0 % UT for 1 cycle 70 % UT for 25/30 cycles (50/60 Hz) 0 % UT for 250/300 cycles (50/60 Hz)	NOT APPLICABLE	Icare HOME (TA022) tonometer is not operational when connected to an external computer, Icare HOME (TA022) batteries are not rechargeable
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.  WARNING: Sources of power frequency magnetic field should be used no closer than 15 cm (6 inches) to any part of Icare HOME (TA022), including cables specified by the manufacturer. Otherwise, degradation of the performance could result.

Guidance and manufacturer's declaration – Electromagnetic immunity			
Icare HOME (TA022) is intended for use and should be used in the electromagnetic environment specified below.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment-Guidance
Conducted RF IEC 61000-4-6	3 V 0,15 MHz – 80 MHz	3 V	 WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Icare HOME (TA022), including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
	6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	6 V	
Radiated RF IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz	10 V/m	
Proximity fields from RF wireless communications equipment IEC 61000-4-3	380 - 390 MHz 27 V/ m; PM 50%; 18 Hz	27 V/m	Interference may occur in the vicinity of equipment marked with the following symbol: 
	430 - 470 MHz 28 V/m; (FM ±5 kHz, 1 kHz sine) PM; 18 Hz	28 V/m	
	704 - 787 MHz 9 V/m; PM 50%; 217 Hz	9 V/m	
	800 - 960 MHz 28 V/m; PM 50%; 18 Hz	28 V/m	
	1700 - 1990 MHz 28 V/m; PM 50%; 217 Hz	28 V/m	
	2400 - 2570 MHz 28 V/m; PM 50%; 217 Hz	28 V/m	
	5100 - 5800 MHz 9 V/ m; PM 50%; 217 Hz	9 V/m	

21. TRAINING

All health care professionals, patients and caregivers should practice the use of the tonometer with this instruction manual or the patient guide and a training video, or receive training from a trained professional prior to performing tonometry.

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